Collection Calendar

Insert the days of week starting with the first day in which your menstrual cycle began. ie. Tues, Weds, Thurs, etc.



Your collection schedule is based on the first day of your menstrual cycle (day menstrual blood flow begins). Use this calendar to log the following dates to assist your healthcare provider with interpretation of your results: Enter onset date of **PREVIOUS** menstrual period enter onset date of **PREVIOUS** menstrual period all collection dates in the chart below; and onset date of **PREVIOUS** menstrual period (onset date of **PREVIOUS** menstrual period (o

*Note: this calendar is meant to serve as a quide to aid in your collection. See full kit instructions for details on acceptable specimen collection.

Day 1	Day 2	Day 3	Day 4	Day 5	Day 6	Day 7
Do not collect First day of your menstrual cycle		1st Collection 7-9 AM Label #1		2nd Collection 7-9 AM Label #2		
Date	Date	Date	Date	Date	Date	Date
Day 8	Day 9	Day 10	Day 11	Day 12	Day 13	Day 14
3rd Collection 7-9 AM Label #3			4th Collection 7-9 AM Label #4	5th Collection 7-9 AM Label #5		6th Collection 7-9 AM Label #6
Date	Date	Date	Date	Date	Date	Date
Day 15	Day 16	Day 17	Day 18	Day 19	Day 20	Day 21
	7th Collection 7-9 AM Label #7		8th Collection 7-9 AM Label #8		9th Collection 7-9 AM Label #9	
Date	Date	Date	Date	Date	Date	Date
Day 22	Day 23	Day 24	Day 25	Day 26	Day 27	Unknown
	10th Collection 7-9 AM Label #10			11th Collection 7-9 AM Label #11		12th Collection First day of your menstrual cycle Label #12
Date	Date	Date	Date	Date	Date	Date



Collection Schedule



DOB

Sample 4

Day 11

This form must be completed and returned with your samples

Your collection schedule is based on the first day of your menstrual cycle (the day menstrual blood flow begins). You may find it helpful to use a calendar to remind you of collection dates.

If you are not experiencing periods, you may start on any day, however if you have a recurring menstrual symptom, try to target the test so that the symptom falls on day 11 or 12. The number of samples will help to identify any underlying patterns.

Collect your saliva during the	Day 1	Count from this day (the first day of your menstrual cycle). Do not collect on this day.	
specified time frame	Day 2	No collection on this day.	
each day:	Days 3, 5, 8, 11, 12, 14, 16, 18, 20, 23, 26	Collect the saliva samples in the morning, and at approximately the same time for each sample.	Use labels #1, 2, 3, 4 etc.
	First day of your NEXT period		Label #12

Please ensure you have read the instructions fully prior to commencing collection of samples.

Key points to remember:

- Avoid alcohol for 12 hours prior to collecting each sample.
- One hour prior to collection do not eat, brush or floss your teeth, use mouthwash, chew gum or use any tobacco products. You may drink ONLY water during this time.
- If you make a mistake or need to restart collecting the samples, please rinse the tube(s) with water only, and allow them to dry naturally.

		'	— ' — ' —				
Sample 5 Day 12 / /	Sample 6 Day 14 / /	Sample 7 Day 16 / /	Sample 8 Day 18 / /				
Sample 9 Day 20//	Sample 10 Day 23//	Sample 11 Day 26//	Sample 12 Day 28				
Questions The following questions MUST be completed and returned with your samples.							
Please answer ALL the following questions by circling the appropriate response: $ \\$							
Do you have a re	cle? Y	ES NO					
What is your average cycle length? (e.g. 28 days) To Days							
• Do you feel you may be entering menopause? YES NO							
Do you experien	PMS? Y	ES NO					
Is this test being	sons? Y	ES NO					
If any hormones or drugs have been/are being taken within 6 months of this test, please indicate below.							
Progesterone (Oral	Last Take	Last Taken					
Oestrone / Oestrad	Last Take	Last Taken					
Testosterone	Last Take	Last Taken					
DHEA	Last Take	Last Taken					
Other (specify)	Last Take	Last Taken					

Surname

Sample 3

Day 8

Note the actual dates of collection below:

Sample 2

Day 5

Patient details: First name

Sample 1

Day 3 (First)



©2017 Genova Diagnostics Europe